



BURETTES & PIPETTES


PIPETTE, MEDICINE DROPPER STRAIGHT GLASS WITH RUBBER BULB

	(70)	BP 170	3"		\$.30
		BP 170-12		dozen	3.50
		BP 172	4"		.35
		BP 172-12		dozen	3.75


PIPETTE, MEDICINE DROPPER BENT TIP 3" WITH RUBBER BULB

	(70)	BP 171			\$.30
		BP 171-12		dozen	3.50


PIPETTE, MEDICINE DROPPER STRAIGHT PLASTIC 3" WITH RUBBER BULB

		BP 173			.30
		BP 173-12		dozen	3.50


DROPPING PIPETTE, GLASS TUBE WITH RED RUBBER BULB

		BP 204	4"		\$ 1.50
		BP 206	6"		1.60
		BP 208	8"		1.70
		BP 210	10"		1.80
		BP 212	12"		1.90

PROTOZOA PIPETTE, LARGE VOLUME, GLASS TUBE WITH BULB

		BP 220	200 mm X 30 mm		\$13.95
---	--	--------	----------------	--	---------

PLASTIC PIPETTE, 5 ml ALL IN ONE OF FLEXIBLE PLASTIC W/ SQUEEZE BULB, GOOD FOR CARTISIAN DIVER EXPERIMENT

		BP 305			\$.15
		BP 305-12		dozen	1.50

**VOLUMETRIC PIPETTE, GLASS WITH ONE GRADUATE, ACCURATE AT 25° C
CALIBRATED TO DELIVER**

BP 401	1 ml	\$4.50
BP 405	5 ml	4.50
BP 406	10 ml	5.25
BP 409	25 ml	7.35
BP 410	50 ml	9.00



PIPETTE, MOHR, GRADUATED OF BOROSILICATE GLASS

(925)


CAPACITY SUBDIVISIONS

BP 421	1 ml	1/100 ml	\$1.20
BP 422	2 ml	1/100 ml	1.35
BP 425	5 ml	1/10 ml	1.45
BP 426	10 ml	1/10 ml	1.95
BP 429	25 ml	2/10 ml	2.19



BURETTE, GRADUATED FOR LIQUIDS OF PYREX GLASS

(29)	BP 509	25 ml	\$ 9.00
	BP 512	50 ml	10.98
	BP 520	100 ml	12.50



**BURETTE ATTACHMENT CONSISTING OF GLASS TIP, PINCH COCK & PIECE OF
RUBBER TUBING**



BP 530 \$ 2.25

BURETTE WITH GLASS STOPCOCK, GRADUATED

(948)	BP 606	10 ml	\$ 15.50
	BP 609	25 ml	18.50
	BP 612	50 ml	22.90
	BP 620	100 ml	25.20



BURETTES & PIPETTES

BURETTE WITH TEFLON STOPCOCK, GRADUATED

(948-A)	BP 659	25 ml	\$32.50
	BP 662	50 ml	34.00
	BP 670	100 ml	36.10



PIPETTE FILLER, CHEMICAL RESISTANT BULB, PRESS TO FILL - PRESS TO RELEASE FITS ANY SIZE PIPETTE

(8610)	BP 840	\$ 8.00
--------	--------	---------




PIPETTE PUMP, FAST RELEASE PROVIDES FOR EASY & PRECISE PIPETTING. W/ QUICK RELEASE LEVER & THUMB WHEEL, SOFT ELASTIC CHUCK W/THREADED COLAR TO TIGHTLY GRIP GLASS OR PLASTIC PIPETTES, COLOR CODED TO SIZE

BP 802	2 ml	\$ 8.50
BP 806	10 ml	9.50
BP 809	25 ml	10.00



PIPETTE SUPPORT, "A" FRAME STYLE WITH HORIZONTAL STORAGE FOR 12 PIPETTES, POLYPROPYLENE

BP 852	\$7.95
--------	--------



PIPETTE SUPPORT, ROUND TYPE WITH VERTICAL STORAGE FOR 28 PIPETTES

BP 854	\$29.95
--------	---------

