

CLAMPS

SCREW CLAMP, HEAVY NICKEL PLATED, CAN BE USED TO STOP OR CONTROL FLOW IN RUBBER OR PLASTIC TUBING, OPENING SIZE 1/2"

(7350)

CL. 115

\$1.25



HOFFMAN TUBING CLAMP, NICKEL PLATED FINISH, OPENING SIZE 3/4"

(7221)

CL. 120

\$1.50



CLAMP, PINCH TYPE (MOHR) TO STOP FLOW IN RUBBER OR PLASTIC TUBING. OF NICKEL PLATED BRASS.

(43)

CL. 125

\$1.30



CLAMP, PINCH (DAYS) OF STAINLESS STEEL.

(43D)

CL. 130

\$1.45




CLAMP, HOFFMAN SCREW, W/SWING BAR FOR EASIER INSERTION OF TUBING, PLATE FINISH, OPENING SIZE 1/2"

(722X)

CL. 135

\$1.30



CLAMP, HOFFMAN, SINGLE SCREW OF BRASS AND NICKEL PLATE CAN BE USED TO SLOW-UP OR STOP FLOW IN TUBING, OPENING SIZE 1/2"

(735)

CL. 140

\$1.25



LAMPS

LAMP, UTILITY 3 FINGER JAWS, MEDIUM VINYLIZED, FOR HOLDING IRREGULAR SHAPED GLASSWARE UP TO 2½" DIAMETER

(734) CL 145

\$9.99



LAMP, TEST TUBE, HEAVILY PLATED SPRING STEEL.

(44) CL 150

\$1.35



LAMP, TEST TUBE (STODDARD) BEST QUALITY, WITH "FINGER GRIPS".

(44B) CL 155

\$1.59



LAMP UTILITY "NEVER RUST" JAWS OPEN TO 2½", MEDIUM SIZE

(2733) CL 205

\$5.95



LAMP, THERMOMETER ADJUSTABLE 360° JAW

CL 220

\$ 5.95



CLAMP, WALL-TYPE WITH SCREW TO HOLD TUBES, BURETS, ETC. TO WALL OR SHELF BOARD

(3-36)

CL 236

\$5.20




CLAMP BURETTE "NEVER RUST" HOLDS FLASKS, BURETTES, OR TUBES.

(45)

CL 250

\$6.75



CLAMP BURETTE, PLASTIC SPRING LOADED

CL 251

CL 252

SINGLE

\$4.50

DOUBLE

6.95



CLAMP BURETTE, LINCOLN, DOUBLE STYLE MADE ENTIRELY OF STAMPED STEEL, PLATED FINISH, HOLDS 2 BURETTES

(729)

CL 255

\$6.95



CLAMP BURETTE, SMALL FOR USE WITH 18" ROD RINGSTANDS ONLY, ALL METAL, PLATED FINISH

(3-71)

CL 271

\$2.75



CLAMPS

CLAMP EXTENSION SMALL, HOLDS ITEMS 1/32" TO 1" DIAMETER "NEVER RUST"
(82 A) CL 300 \$4.50



CLAMP EXTENSION, OF PLATED METAL AND COATED JAWS, JAWS OPEN
TO 2 1/4"

(45 A) CL 305 \$5.75



CLAMP EXTENSION UNIVERSAL, TRIPRONG JAWS, COATED FOR HOLDING
LARGER GLASSWARE, CONDENSERS, ETC., JAWS OPEN TO 3 1/4".

(45 B / 2726) CL 307 \$6.89



CLAMP HOLDER, TO HOLD ABOVE EXTENSION CLAMPS AND RINGS.

(45 C) CL 310 REGULAR \$4.55



CL 312 LARGE HEAVY DUTY \$4.99



(4746) CL 311 ADJUSTABLE \$5.50



CLAMP HOLDER, CONTORT HOLDS CLAMPS OR RINGS AT ALMOST ANY ANGLE
OF "NEVER RUST" ALLOY

(2750) CL 315 \$10.50

