

FUNNELS, THISTLE TUBES & FILTER PAPER

FILTER PAPER FOR BUCHNER AND OTHER FUNNELS, MEDIUM FLOW (PK OF 50)



(59)

FP 155	55 mm	\$ 1.25
FP 170	70 mm	1.50
FP 109	90 mm	1.95
FP 111	110 mm	2.20
FP 112	125 mm	2.75
FP 115	150 mm	3.50
FP 118	180 mm	4.25

FUNNEL, BUCHNER POLYPROPYLENE, 2 PIECE SEPARATES FOR EASY CLEANING
(253 P)




FILTER PAPER SIZE

FU 251	55 mm	\$ 7.30
FU 253	70 mm	9.00
FU 255	90 mm	11.85
FU 256	110 mm	19.50
FU 257	150 mm	25.80

FUNNEL, BUCHNER PORCELAIN 1 PIECE CONSTRUCTION

(235B)



FU 261	55 mm	\$6.75
FU 263	70 mm	7.75
FU 265	90 mm	8.90
FU 266	110 mm	10.20
FU 267	125 mm	14.50

FUNNEL, SHORT STEM, HEAVY BOROSILICATE GLASS, STEM LENGTH SAME AS DIAMETER



FU 424	40 mm	\$2.60
FU 425	50 mm	3.50
FU 427	75 mm	4.10
FU 429	90 mm	4.50
FU 430	120 mm	5.20

FUNNELS, THISTLE TUBES & FILTER PAPER

FUNNEL TUBE, CONE TOP, 150 mm STEM

(63C)

FU 510

\$ 2.75



FUNNEL, BURETTE FILLER, CHONICAL TOP 40 mm WITH 50 mm STEM

(20)

FU 520

\$ 1.59



FUNNEL TUBE, THISTLE TOP, DESIGNED TO PREVENT BACK SPLASH WHEN
POURING CORROSIVE MATERIALS

(63P)

FU 525

150 mm STEM

\$6.00

FU 526

300 mm STEM

6.00

FU 527

400 mm STEM

6.00



FUNNEL, THISTLE TUBE WITH BULB TRAP OF PYREX

(3-21)

FU 530

\$7.95



FUNNEL GLASS, 40 mm (1.5")

(62A)

FU 550

\$3.99



**FUNNEL, CHEMICAL OF MOLD BLOWN GLASS NOMINAL DIAMETER 55 mm
STEM LENGTH 65 mm**

(62)

FU 552

\$3.99



FUNNEL, BUNSEN OF CLEAR GLASS, STEM LENGTH 150 mm

(90)

FU 554

100 mm

\$ 5.95

FU 555

90 mm

4.95

FU 556

65 mm

4.00

FU 557

55 mm

3.80

FU 558

45 mm

3.20



FUNNEL, CONVOLUTED INSIDE FOR EASY FILTER FLOW , OF CLEAR GLASS

FU 560

35 mm

\$ 3.50

FU 561

45 mm

3.75

FU 562

55 mm

4.50

FU 563

100 mm

6.50



**FUNNEL CHEMICAL HAND BLOWN OF HEAVY CLEAR LAB GLASS FOR POURING
AND FILTERING 16 oz CAPACITY**

(5700)

FU 570

\$6.50



FUNNEL, THISTLE TUBE, WITH SQUARE MOLDED TOP TO PREVENT ROLLING.

(5702)

FU 572

350 mm STEM

\$4.50



FUNNELS, THISTLE TUBES & FILTER PAPER

FUNNEL, THISTLE TUBE NALGENE PLASTIC 300 mm STEM

(942) FU 580

\$3.50



FUNNEL, POWDER WITH SHORT WIDE STEMS, PYREX GLASS

TOP DIAMETER



FU 590	65 mm	\$ 8.99
FU 591	75 mm	9.99
FU 592	100 mm	14.69
FU 593	125 mm	29.90
FU 594	150 mm	57.99

FUNNEL, POLYPROPYLENE PLASTIC, RIBBED AND VENTED 65 mm DIAMETER

(62P) FU 620

\$1.35



FUNNEL, POLYPROPYLENE, 3 EXTERNAL RIBS, STEM LENGTH APPROXIMATELY SAME AS DIAMETER.



FU 635	55 mm	\$.65
FU 636	65 mm	.85
FU 637	75 mm	1.10
FU 639	100 mm	1.85
FU 640	150 mm	2.45

FUNNEL, POLYPROPYLENE, 8 INTERNAL RIBS TO AID IN FILTRATION, 32 mm STEM LENGTH.



FU 683	35 mm	\$.39
FU 685	52 mm	.49
FU 686	65 mm	.59

FUNNEL, SHORT STEM**POLYETHYLENE**

FU 696 65 mm

\$1.20

FU 698 125 mm

1.50

POLYPROPYLENE

FU 699 85 mm

\$1.35

FUNNEL, SEPARATORY GLOBE SHAPE OF PYREX GLASS, GLASS STOPPERED

(9464)



FU 720 50 ml \$ 26.95

FU 722 100 ml 31.95

FU 724 250 ml 37.95

FU 727 500 ml 39.50

FU 728 1000 ml 49.90

FU 729 2000 ml 69.90

FUNNEL SEPARATORY SQUIBB, PYREX WITH GLASS STOPPER AND STOPCOCK

(943)



FU 751 100 ml \$ 19.95

FU 752 125 ml 22.95

FU 754 250 ml 25.95

FU 757 500 ml 32.95

FU 760 1000 ml 54.95

FU 762 2000 ml 155.50

RUBBER FILTER AID, FLEXIBLE RUBBER PLATE TO SIT ON TOP OF FLASK AND SEAL FUNNEL TO FLASK. HOLDS MOST SIZES OF FUNNELS

FU 910

\$ 4.95

