

## GLASSWARE WITH GROUND GLASS JOINTS

24/40 JOINTS UNLESS NOTED

OTHER SIZES & ITEMS ARE AVAILABLE, REQUEST A QUOTE

### FLASK ERLLENMEYER, 24/40 JOINT



| (794) | GC 401 | 125 ml  | \$ 23.00 |
|-------|--------|---------|----------|
|       | GC 402 | 250 ml  | 25.00    |
|       | GC 405 | 500 ml  | 29.00    |
|       | GC 410 | 1000 ml | 35.00    |
|       | GC 406 | 2000 ml | 45.00    |
|       | GC 407 | 4000 ml | 79.00    |
|       | GC 408 | 6000 ml | 110.00   |

### FLASK ROUND BOTTOM, 24/40 JOINT



| (945) | GC 425 | 50 ml   | \$19.00 |
|-------|--------|---------|---------|
|       | GC 426 | 100 ml  | 20.00   |
|       | GC 424 | 250 ml  | 22.50   |
|       | GC 427 | 500 ml  | 29.95   |
|       | GC 430 | 1000 ml | 39.95   |
|       | GC 440 | 2000 ml | 49.95   |
|       | GC 441 | 3000 ml | 65.00   |
|       | GC 439 | 5000 ml | 98.00   |

### FLASK TWO NECK ROUND BOTTOM, 24/40 JOINTS



| (142) | GC 444                                  | 250 ml  | \$ 50.00 |
|-------|---|---------|----------|
|       | GC 445 <td>500 ml</td> <td>52.00</td>   | 500 ml  | 52.00    |
|       | GC 442 <td>1000 ml</td> <td>55.50</td>  | 1000 ml | 55.50    |
|       | GC 446 <td>2000 ml</td> <td>75.00</td>  | 2000 ml | 75.00    |
|       | GC 443 <td>3000 ml</td> <td>85.00</td>  | 3000 ml | 85.00    |
|       | GC 447 <td>5000 ml</td> <td>120.00</td> | 5000 ml | 120.00   |

### FLASK TWO NECK, ANGLED SIDE NECK AT AN ANGLE TO CENTER 24/40 JOINTS



|  |        |         |          |
|--|--------|---------|----------|
|  | GC 380 | 250 ml  | \$ 60.00 |
|  | GC 381 | 500 ml  | 65.00    |
|  | GC 382 | 1000 ml | 75.00    |
|  | GC 383 | 2000 ml | 90.00    |
|  | GC 384 | 3000 ml | 100.00   |

### FLASK THREE NECK, 24/40 JOINTS




| (801) | GC 412                                  | 250 ml  | 58.00  |
|-------|---|---------|--------|
|       | GC 415 <td>500 ml</td> <td>65.00</td>   | 500 ml  | 65.00  |
|       | GC 417 <td>1000 ml</td> <td>79.00</td>  | 1000 ml | 79.00  |
|       | GC 420 <td>2000 ml</td> <td>119.00</td> | 2000 ml | 119.00 |
|       | GC 423 <td>3000 ml</td> <td>130.00</td> | 3000 ml | 130.00 |
|       | GC 421 <td>5000 ml</td> <td>170.00</td> | 5000 ml | 170.00 |

### FLASK THREE NECK, ANGLED SIDE NECKS AT AN ANGLE TO CENTER 24/40 JOINTS



|  |        |         |          |
|--|--------|---------|----------|
|  | GC 390 | 250 ml  | \$ 68.00 |
|  | GC 391 | 500 ml  | 75.00    |
|  | GC 392 | 1000 ml | 85.00    |
|  | GC 393 | 2000 ml | 125.00   |
|  | GC 394 | 3000 ml | 135.00   |

**CONDENSER ALLIHN, BOTH ENDS 24/40 JOINTS**



|       |        | <b>JACKET LENGTH</b> |  |         |
|-------|--------|----------------------|--|---------|
| (410) | GC 480 | 200 mm               |  | \$60.00 |
|       | GC 481 | 250 mm               |  | 65.00   |
|       | GC 482 | 300 mm               |  | 69.00   |
|       | GC 483 | 400 mm               |  | 75.00   |
|       | GC 484 | 500 mm               |  | 80.00   |
|       | GC 485 | 600 mm               |  | 85.00   |
|       | GC 486 | 750 mm               |  | 92.00   |
|       | GC 487 | 900 mm               |  | 110.00  |
|       | GC 488 | 1200 mm              |  | 210.00  |


**CONDENSER GRAHAM, BOTH ENDS 24/40 JOINTS**



|       |        | <b>JACKET LENGTH</b> |  |          |
|-------|--------|----------------------|--|----------|
| (460) | GC 490 | 200 mm               |  | \$ 80.00 |
|       | GC 491 | 250 mm               |  | 80.00    |
|       | GC 492 | 300 mm               |  | 85.00    |
|       | GC 493 | 400 mm               |  | 98.00    |
|       | GC 494 | 500 mm               |  | 110.00   |
|       | GC 495 | 600 mm               |  | 125.00   |
|       | GC 496 | 700 mm               |  | 145.00   |
|       | GC 497 | 900 mm               |  | 165.00   |
|       | GC 498 | 1200 mm              |  | 200.00   |

**CONDENSER WEST**


(550)



|        |  |          |
|--------|--|----------|
| GC 454 | 24/40 JOINT ON DRIP END, 300 mm JACKET   | \$ 37.75 |
| GC 456 | 24/40 JOINTS ON BOTH ENDS, 300 mm JACKET | 55.00    |
| GC 458 | 24/40 JOINTS ON BOTH ENDS, 400 mm JACKET | 59.75    |

**CONDENSER, FRIEDRICHS 24/40 JOINTS, SPIRAL INNER TUBE, WITH INNER T/S JOINT AT BOTTOM AND OUTER T/S JOINT IN SIDE ARM**

(9419) GC 460 \$100.00



**CONDENSER DRIP TUBE, 24/40 JOINT**

(265)

GC 462

\$25.00



## GLASSWARE WITH GROUND GLASS JOINTS

EXTRACTION APPARATUS, SOXHLET WITH GROUND GLASS 24/40 JOINTS. CONSISTS OF ONE EACH: CONDENSER, EXTRACTION TUBE AND FLASK. THE FLASK CAPACITIES ARE MEASURED TO THE BASE OF THE NECK. COMPLETE AS ILLUSTRATED



(5158)

|        | FLASK SIZE |          |
|--------|------------|----------|
| GC 470 | 125 ml     | \$235.00 |
| GC 471 | 250 ml     | 245.00   |
| GC 472 | 500 ml     | 285.00   |

## SEPARATORY FUNNEL, 24/40 JOINTS



(411)

|        |         |          |
|--------|---------|----------|
| GC 521 | 125 ml  | \$ 65.00 |
| GC 524 | 250 ml  | 79.00    |
| GC 525 | 500 ml  | 99.00    |
| GC 527 | 1000 ml | 170.50   |
| GC 529 | 2000 ml | 219.00   |

## DISTILLING COLUMN VIGREUX WITH INDENTATIONS OVER THE ENTIRE LENGTH FOR EXCELLENT VAPOR-LIQUID CONTACT



|        |        |          |
|--------|--------|----------|
| GC 530 | 150 mm | \$ 56.00 |
| GC 531 | 200 mm | 58.00    |
| GC 533 | 300 mm | 73.00    |
| GC 535 | 450 mm | 106.00   |
| GC 537 | 600 mm | 146.00   |

**SAFETY CLIPS FOR STANDARD 24/40 JOINTS. HOLD JOINTS SECURLY TOGETHER,  
SAFE FOR TEMPERATURES TO 150° C, EXCELLENT CHEMICAL RESISTANCE**

GC 540

\$ 3.50



**STOPPER, GLASS PENNYHEAD 24.40 JOINT**

GC 545

\$ 16.50



**ADAPTER CLAISEN 24/40 JOINTS**

(161)

GC 550


\$69.00



**ADAPTER CLAISEN W/ THERMOMETER FITTING 24/40 JOINTS**

GC 565

\$ 59.00



**ADAPTER, CLAISEN W/ # 7 THREADED TOP FOR USE WITH A THERMOMETER**

GC 567

\$ 46.00



**ADAPTER, CLAISEN W/ TAKE OFF STOPCOCK 2 mm STOPCOCK WITH TEFLON PLUG**

GC 568


\$ 69.00



**ADAPTER, CLAISEN MODIFIED TOP OF SIDE TUBE HAS A #7 THREAD WITH BUSHING  
TAKE OFF TUBE HAS AN INNER DRIP TUBE**

GC 569

\$ 59.00



## GLASSWARE WITH GROUND GLASS JOINTS

ADAPTER, CONNECTING TUBE, 75° ANGLE, MALE 24/40 JOINTS



(409)

GC 548

\$ 25.50

ADAPTER, CONNECTING TUBE 3-WAY 75° ANGLE, ALL 24/40 JOINTS



(851)

GC 554


\$42.50



GC 555

45.00


ADAPTER, CONNECTING TUBE 105° ANGLE, MALE & FEMALE 24/40 JOINTS



GC 572

\$ 45.00

ADAPTER, CONNECTING TUBE 200 mm CENTER TO CENTER AT 75° INNER JOINTS  
24/40



GC 574

\$ 37.00

ADAPTER, CONNECTING TUBE 90° BEND 24/40



GC 576

\$ 26.00

ADAPTER, CONNECTING TUBE OFFSET 24/40



GC 578

\$ 25.00

ADAPTER, CONNECTING TUBE STRAIGHT 24/40



GC 580

30 mm

\$ 23.00

GC 582

70 mm

24.00

GC 584

120 mm

26.00

GC 586

175 mm

27.00

ADAPTER, CONNECTING TUBE INNER JOINTS 24/40 STRAIGHT



GC 590

30 mm

\$ 25.00

GC 591

70 mm

27.00

GC 592

120 mm

29.00

GC 593

170 mm

32.00

**ADAPTER, CONNECTING TUBE OUTER JOINTS 24/40, STRAIGHT**

|        |        |          |
|--------|--------|----------|
| GC 594 | 30 mm  | \$ 26.00 |
| GC 595 | 70 mm  | 28.00    |
| GC 596 | 120 mm | 30.00    |
| GC 597 | 170 mm | 33.00    |

**ADAPTER, DISTILLATION 75° WITH INTERNAL #7 THREAD AND BUSHING FOR THERMOMETER OR BLEED TUBE**

|        |          |
|--------|----------|
| GC 610 | \$ 46.00 |
|--------|----------|

**ADAPTER, DISTILLING TRAP USED TO AVOID CARRY OVER OF REACTIONS WHICH HAVE A TENDENCY TO FOAM**

|        |          |
|--------|----------|
| GC 612 | \$ 41.00 |
|--------|----------|

**ADAPTER, DISTILLATION 105° BEND**

|        |          |
|--------|----------|
| GC 614 | \$ 30.00 |
|--------|----------|

**ADAPTER, DISTILLATION W/ SEALED DRIP TUBE AND HOSE CONNECTION FOR EVACUATION**

|        |          |
|--------|----------|
| GC 616 | \$ 53.00 |
|--------|----------|

**ADAPTER, HOSE END, STRAIGHT 24/40 INNER JOINT, BOTTOM 10mm HOSE TOP**

|              |          |
|--------------|----------|
| (811) GC 552 | \$ 12.98 |
|--------------|----------|

**ADAPTER, HOSE END, 90°**

|        |          |
|--------|----------|
| GC 553 | \$ 13.98 |
|--------|----------|

**ADAPTER, HOSE END W/ GLASS STOPCOCK STRAIGHT**

|        |          |
|--------|----------|
| GC 615 | \$ 45.00 |
|--------|----------|

## GLASSWARE WITH GROUND GLASS JOINTS

ADAPTER, HOSE END WITH GLASS STOPCOCK 90° 24/40

GC 617

\$ 45.00



ADAPTER, INLET THERMOMETER W/ RUBBER ADAPTER FOR BLEED TUBE OR THERMOMETER

GC 626

\$ 24.00



ADAPTER, INLET THERMOMETER W/ # 7 INTERNAL THREADS SUPPLIED WITH BUSHING AND O-RING

GC 627

\$ 30.00



ADAPTER, INLET THERMOMETER / VACUUM 24/40

GC 628

\$ 47.00



ADAPTER, VACUUM STRAIGHT WITH DRIP TUBE 24/40

(420)

GC 556

\$ 62.00



ADAPTER, VACUUM 105° WITH DRIP TUBE

GC 557

\$ 65.00



VACUUM GREASE, TUBE 5.3 OUNCE

(635)

GC 560

\$ 21.00

