

INSTRUMENTS

CORK BORER SET WITH CLEARING PUNCH USED TO CUT HOLES IN RUBBER STOPPERS & CORKS

(50)	IN 150	SET OF 3	3/16", 1/4", 5/16"	\$ 5.00
	IN 152	SET OF 6	4 mm - 12 mm	12.00
	IN 154	SET OF 12	4 mm - 20 mm	24.95



CORK BORER SHARPENER, BRASS WITH WOOD HANDLE AND QUALITY STEEL KNIFE

(946) IN 167 \$6.99



INOCULATING NEEDLE, STRAIGHT METAL HANDLE, NICHROME NEEDLE

(1032) IN 332 \$1.59



INOCULATING LOUPE, METAL HANDLE, NICHROME LOUPE

(1032 L) IN 334 \$1.99



TONGS BEAKER, UNCOATED FOR BEAKERS 50 TO 2000 ml CAPACITY. FLAT RIVETED JOINT. MAXIMUM JAW OPENING 4 1/2"

IN 380 \$ 4.59



TONGS BEAKER, COATED SAME AS ABOVE WITH JAWS COATED WITH PLASTISOL

IN 382 \$ 4.95



TONGS CRUCIBLE OF HEAVY PLATED STEEL, SURE GRIP JAWS

(85) IN 385

\$3.55

**SPATULA WITH WOODEN HANDLE, FLEXIBLE STAINLESS STEEL BLADE WITH ROUNDED END. MOUNTED IN A POLISHED WOOD HANDLE. BY BLADE LENGTH**

IN 413	3"	\$4.20
IN 414	4"	4.50
IN 416	6"	6.50
IN 418	8"	9.00

SPATULA / SPOON STAINLESS STEEL, WILL REACH INTO JARS TO LOOSEN AND SPOON CHEMICALS.

(58)	IN 426	6"	\$2.49
	IN 428	8"	2.99

**SPATULA SEMI MICRO**

(1100) IN 430 \$ 3.00

**SCOOP LABORATORY MADE OF STAINLESS STEEL COMPLETE WITH A REMOVABLE WOOD HANDLE (BY TOTAL LENGTH)**

(1576) IN 440 11" \$2.49



INSTRUMENTS

SPOON, DEFLAGRATION, 3/4" BOWL, 15" HANDLE, STAINLESS STEEL
(76) IN 476 \$3.99




SPOON SEMI MICRO
(2000) IN 478 \$ 3.00



CALIPER VERNIER STEEL DOUBLE SCALE 0 - 120 mm X 0.1 mm AND 5" X 0.25"
ADJUSTMENT ROLLER AND LOCKING SCREW
IN 500 \$ 5.95



CALIPER VERNIER PLASTIC DOUBLE SCALE RANGE 6" OR 150 mm
READABILITY 1/128" OR 0.05 mm
IN 505 \$ 5.95



MICROMETER, SCREW FORM WITH GAUGE 25 X 1 mm PITCH
IN 510 \$ 15.95




SPLINTS, WOOD FOR COMBUSTION IN OXYGEN EXPERIMENTS

(75)	IN 601	BUNDLE OF 25	\$.39
	IN 601 - 50	BUNDLE OF 50	.69
	IN 601 - 500	BUNDLE OF 500	3.50

**STENCIL, OF PLASTIC FOR MAKING DRAWINGS OF APPARATUS SET-UPS**

(78 S)	IN 618		\$1.95
--------	--------	--	--------

**PENCILS "GLASS MARKING"**

(88)	IN 620	RED	\$1.00
	IN 621	BLUE	1.00
	IN 622	BLACK	1.00

**PNEUMATIC TROUGH NALGENE PLASTIC COMPLETE WITH SHELF & OVER-FLOW**

(694)	IN 694	12" X 7" X 3"	\$16.49
-------	--------	---------------	---------

**PNEUMATIC TROUGH GALVANIZED IRON, COMPLETE WITH SHELF & OVER-FLOW**

(694)	IN 695	10" X 7" X 4"	\$21.95
-------	--------	---------------	---------

