

PHYSICAL SCIENCE

LABORATORY JACK 15 X 12 cm TOP PLATE & BASE, VERTICAL HEIGHT 26 cm

ALSO AVAILABLE WITH REMOVABLE 20 X 20 cm PLATFORM & 45 cm SUPPORT ROD

PY 440	15 X 12 cm PLATFORM	\$59.00
PY 444	20 X 20 cm PLATFORM	65.00



LENSES, UNMOUNTED, SIZE AND FOCAL LENGTH IN mm, EACH


DIAMETER	FOCAL LENGTH	DOUBLE CONVEX	DOUBLE CONCAVE
38	50	PY 500 \$ 2.00	—
38	100	PY 502 1.50	PY 527 \$ 2.00
38	150	PY 503 1.50	PY 528 2.00
38	200	PY 504 1.50	PY 529 2.00
38	300	PY 505 1.50	—
50	50	PY 506 3.00	PY 531 3.50
50	100	PY 508 2.50	PY 533 2.75
50	150	PY 507 2.50	PY 532 2.75
50	200	PY 509 2.50	PY 534 2.75
50	250	PY 510 2.50	—
75	175	PY 515 3.50	—
75	200	—	PY 540 4.50
100	200	PY 519 6.00	—
100	500	PY 521 6.00	—



LENS SET, DEMONSTRATION 6 LENS SET, 1 EACH: DOUBLE CONVEX, PLANO CONVEX, CONCAVO-CONVEX, DOUBLE CONCAVE, PLANO CONCAVE, CONVEXO-CONCAVE, 75 mm DIAMETER

PY 550

\$11.25



PHYSICAL SCIENCE

POCKET LOUPE 10X MAGNIFICATION, COMPACT FOLDING LOUPE, ALL METAL CASE
 PY 560

\$5.75



MAGNIFIER FOLDING WITH POUCH 50 mm DIAMETER 2.5 X MAGNIFICATION
 PY 561

\$1.99



MAGNIFIER FOLDING PLASTIC CASE

PY 562	SINGLE LENS	5 X MAG	\$2.00
PY 563	DOUBLE LENS	10 X MAG	3.00
PY 564	TRIPLE LENS	15 X MAG	3.50



MAGNIFIER TRIPOD, 10X, SCREW FOCUS, METAL BASE
 (915) PY 567

\$5.99




MAGNIFIER HAND, 3X, METAL RIMS AND HANDLES, GLASS LENS

PY 570	50 mm	\$ 3.65
PY 571	65 mm	3.95
PY 572	75 mm	4.75
PY 573	100 mm	5.65



PULLEYS 50 mm DIAMETER PLASTIC SHEAVES WITH DEEP V-GROOVES, METAL SWIVEL HOOKS ON BOTH ENDS

PY 601	SINGLE	\$1.25
PY 602	DOUBLE	2.00
PY 603	TRIPLE	2.75
PY 604	QUAD	3.00



PULLEY WITH CLAMP SINGLE 50 mm DIAMETER PULLEY MOUNTED ON A CLAMP TO ATTACH TO SURFACES UP TO 25 mm THICK

PY 605		\$4.20
--------	--	--------



PRISM, EQUILATERAL POLISHED OF FINEST FLINT GLASS

(947)	PY 705	50 mm	\$ 5.25
	PY 707	75 mm	5.99




RADIOMETER, SOLAR POWERED ENGINE

(445)	PY 725		\$ 8.95
-------	--------	--	---------



TUNNING FORKS MADE OF ALUMINUM ALLOY, FREQUENCY & SCALE STAMPED ON EACH FORK, FREQUENCY AVAILABLE: CPS 100 128 256 288 320 341.3 384 426.6 480 512 1024 2048 4096

PY 751	SPECIFY FREQUENCY	\$ 3.90
PY 750-13	SET OF 13	39.90



PHYSICAL SCIENCE


PUMP VACUUM , HAND OPERATED WITH OR WITHOUT GAUGE THAT READS IN BOTH cm & INCHES OF MERCURY. 15 ml PER STROKE PUMPING RATE, FITTED WITH VACUUM RELEASE VALVE

PY 901	WITH GAUGE	\$26.00
PY 902	WITHOUT GAUGE	16.50



PUMP PRESSURE / VACUUM, VACUUM TO 20" Hg, PRESSURE TO 18 psig. CORROSION RESISTANT, NONCONTAMINATING PLASTIC HEAD

PY 904		\$179.00
--------	--	----------



PUMP FILTER, ASPIRATOR, METAL, CONNECTS TO A FAUCET TO PROVIDE VACUUM FOR FAST FILTERING, DISTILLING. COMPLETE WITH FAUCET FITTING, 3/4" NPT

(S9 A)	PY 910	\$15.99
--------	--------	---------

