

STIRRERS, HOT PLATES & BATHS

MAGNETIC STIRRER COMPACT, HI 190M, 1 LITER CAPACITY, SPEED RANGE 100-1000 RPM, 110V POWER, ABS PLASTIC CASE 4.8"X4.8"X1.8" 1.4 LB

SH 190

\$ 76.50



MAGNETIC STIRRER W/STAINLESS STEEL COVER, COMPACT, HI 200M, 1 LITER CAPACITY, SPEED RANGE 100-1000 RPM, 110V POWER, ABS PLASTIC CASE 4.8"X4.8"X1.8" 1.4 lb

SH 200

DISCONTINUED



MAGNETIC STIRRER HEAVY DUTY, HI 300N, 2.5 LITER CAPACITY, SPEED RANGE 100-1000 RPM, 110V POWER, ABS PLASTIC CASE WITH STAINLESS STEEL COVER, 7.1"X7.1"X2.8" 3.6 lb

SH 300

DISCONTINUED



HOT PLATE FOR RINGSTAND, USE SUPPORT PLATE AS A HOT PLATE .
 CONTAINERS ARE PLACED ON A 12 cm (4½") DIAMETER CAST-IRON PLATE WHICH
 CAN BE HEATED WITH A BUNSEN BURNER. A 4.8 mm (3/16") RAISED RIM KEEPS
 CONTAINERS FROM SLIPPING OFF ACCIDENTALLY. PLATE CAN BE LOCKED
 FIRMLY TO SUPPORT RODS OR LATTICES WITH INTEGRAL CLAMP


(5400) SH 400

\$ 11.99



HOT PLATE, CERAMIC TOP IDEAL FOR LOW TO MODERATE TEMPERATURE
 APPLICATIONS. RANGE OF 5° TO 550° C (41° TO 1022° F) FAST HEAT WILL REACH
 MAXIMUM TEMPERATURE WITHIN 10 MINUTES. TEMPERATURE IS INFINITELY
 ADJUSTABLE WITH A PERCENTAGE INPUT CONTROLLER. HEAT INDICATOR
 LIGHT SHOWS WHEN POWER IS APPLIED. DURABLE CERAMIC TOP REMAINS FLAT
 EVEN AT HIGH TEMPERATURES AND IS EASY TO CLEAN. LOW PROFILE DESIGN
 MAKES UNIT STABLE AND EASY TO USE. THREE SIZES AVAILABLE.

	TOP PLATE SIZE	HEATING AREA	
SH 411	4" x 4"	16 SQ IN	\$ 225.00
SH 412	7" x 7"	49 SQ IN	250.00
SH 413	10" x 10"	108 SQ IN	395.00



HOT PLATE / STIRRER CERAMIC TOP STIRRING SPEEDS FROM 60 TO 1200 rpm.
 QUIET RUNNING DIRECT DRIVE MAGNETIC COUPELED ASSEMBLY PROVIDES
 ENOUGH TORQUE TO PULL A VORTEX IN TWO LITERS OF WATER. MAXIMUM
 TEMPERATURE REACHED IN 10 MINUTES. TEMPERATURE RANGE 5° TO 550° C
 INFINITELY ADJUSTABLE. DURABLE CERAMIC TOP REMAINS FLAT EVEN AT HIGH
 TEMPERATURES AND IS EASY TO CLEAN. OVERALL LOW PROFILE MAKES UNIT
 STABLE AND EASY TO USE. THREE SIZES AVAILABLE.

	TOP PLATE SIZE	HEATING AREA	
SH 421	4" x 4"	16 SQ IN	\$ 390.00
SH 422	7" x 7"	49 SQ IN	418.00
SH 423	10" x 10"	108 SQ IN	525.00



WATER BATH IDEAL FOR HEATING FLASKS, BEAKERS, TUBES, ETC. 8" DIAMETER
 X 2½". CONCENTRIC SIZE RINGS, FLAT BOTTOM, & STEAM ESCAPE OUTLET

SH 602	STAINLESS STEEL	\$79.95
SH 603	ALUMINUM	49.95



STIRRERS, HOT PLATES & BATHS

STIRRING BAR RETRIEVER, MAGNETIC 12" POLY, TO REMOVE STIRRING BAR FROM LIQUIDS



SH 502

\$7.00

STIRRING BAR OCTAGONAL, BAR TYPE, FOR USE WITH MAGNETIC STIRRERS



	DIAMETER	LENGTH	
SH 510	1/8"	1/2"	\$ 4.80
SH 511	5/16"	1/2"	4.80
SH 512	5/16"	5/8"	4.80
SH 513	5/16"	1"	5.10
SH 514	5/16"	1 1/2"	5.40
SH 515	3/8"	1"	15.80
SH 516	3/8"	1 1/2"	15.30
SH 517	3/8"	2"	8.40
SH 518	3/8"	2 1/2"	8.40
SH 519	1/2"	1 1/2"	19.40
SH 520	1/2"	3"	24.00

STIRRING BAR STAR HEAD TYPE, FOR USE WITH MAGNETIC STIRRERS DEEP MOLDED CROSSES ACT AS IMPELLERS PROVIDING A STRONG MIXING VORTEX



	DIAMETER	
SH 525	3/8"	\$ 8.10
SH 526	3/4"	8.80
SH 527	1"	8.85
SH 528	1 1/4"	10.65
SH 529	1 1/2"	13.99

STIRRING BAR EGG SHAPE FOR USE WITH MAGNETIC STIRRERS. OVAL SHAPE REDUCES FRICTION AND INCREASES EFFECTIVENESS IN ROUND BOTTOM FLASKS



	DIAMETER	LENGTH	
SH 535	5/8"	1 1/4"	\$ 10.25
SH 536	3/4"	1 5/8"	14.25
SH 537	3/4"	2 1/2"	16.80

STIRRING BAR DISPOSABLE BAR TYPE FOR USE WITH MAGNETIC STIRRERS



	DIAMETER	LENGTH	
SH 540	3/16"	1/2"	\$ 1.20
SH 541	1/4"	1"	1.80
SH 542	1/4"	1 3/8"	2.10
SH 543	5/16"	2"	3.00