

TEST PAPERS & pH METERS

TEST PAPERS pH PAPER UNIVERSAL FOR TESTING DEGREE OF ACIDITY AND ALKALINITY OF SOLUTIONS, pH RANGE IS 1 TO 14, ACCURACY ± 1.0 Ph
100 STRIPS PER VIAL

TP 100 \$ 1.29

TEST PAPERS VIALS OF 100 STRIPS



TP 109	ACID/NEUTRAL/BASE	\$.65
TP 110	RED LITMUS, pH 4/5-8.3 BASE	.65
TP 111	NEUTRAL LITMUS	.65
TP 112	BLUE LITMUS, pH 4.5-8.3 ACID	.65
TP 113	COBALT CHLORIDE	.65
TP 114	CONGO RED pH 3.0-5.0	.65
TP 115	CRYSTAL VIOLET pH 0.0 - 1.8	.65
TP 116	METHYL ORANGE pH 3.2- 4.4	.65
TP 117	NEUTRAL RED pH 6.8 - 8.0	.65
TP 118	LEAD ACETATE	.65
TP 119	PHENOLPHTHALEIN pH 8.2-10.0	.65
TP 120	POTASSIUM IODIDE, STARCH	.65
TP 121	TURMERIC	.65

TEST PAPER PTC FOR USE WITH HEREDITARY TRAIT STUDIES. A TASTE TRAIT VIAL OF 100

TP 125 \$ 1.00

pH METER: CHECKER 1 RANGE 0.00 TO 14.00 Ph, RESOLUTION 0.1pH, ACCURACY ± 0.2 pH, REPLACEABLE ELECTRODE, BATTERY LIFE 3000 HOURS
CALIBRATION: MANUAL 2 POINT, LCD SCREEN



TP 901		\$34.50
TP 902	SPARE ELECTRODE	18.50

pH METER PHEP 4 RANGE 0.00 TO 14.00pH, 0.0° TO 60.0°C, RESOLUTION 0.1pH, ACCURACY ± 0.1 pH, AUTOMATIC TEMPERATURE COMPENSATION, EXTRACTABLE JUNCTION, WATERPROOF, CARTRIDGE REPLACEMENT ELECTRODES, STABILITY INDICATOR, LARGE LCD SCREEN



TP 911		\$79.00
TP 912	SPARE ELECTRODE	26.00

pH TURTLE TURN YOUR PC INTO A pH METER. PLUGS DIRECTLY INTO SERIAL PORT AND WITH SUPPLIED SOFTWARE BECOMES A RECORDING pH METER RANGE 0.0 TO 14.0 pH, RESOLUTION 0.1pH, ACCURACY ± 0.2 pH, AUTOMATIC TEMPERATURE COMPENSATION

TP 921 \$85.00

pH STORING SOLUTION FOR ELECTRODES

TP 940	1 oz	\$ 1.95
TP 941	4 oz	4.95